



SWISS MADE



## MF 3000-10 SLS

Max tube die diameter

12 mm

0.500"

### Use

Domain of use

Precision medical tubing

Calibrated tubes

Pneumatic tubes

Technical tubes

Compound

PVC / PE / PBT / PA / LSZH / PUR / TPU


# mf 3000

## Characteristics

## MF 3000-IO SLS

MF	Fixed centering
SL	Single Layer
SLS	Single Layer + Skin or Stripe(s)
Distribution mono-flow	Conical distribution cartridge Easy to clean-no holes or slots-no dead spots Constant and smooth compound flow No stress and no shear in the melted compound
Gum Space	Axial die adjustment

## Design



## Options

Connections	Connecting flanges for any type of extruder
By pass	Manual - Pneumatic
Distributor	Special distribution cartridges for any compound
Tooling	Hardened steel - Tungsten carbide Special tooling Design recommendation on request
Accessories	Vacuum adapter Assembly and disassembly cleaning kit
Spare parts	Permanent stock, shipment within 24 hours
Services	Microdia answers your request within 24 hours Technical assistance - Engineering - Commissioning